

Ingredients:

100 g APESIN AP 3 contain 5 g didecyldimethyl ammonium chloride, 6 g cocospropylene diamine guanidinium diacetate, 3.5 g isopropanol, less than 5% non-ionic surfactants. Other ingredients: perfume oil. pH approx. 4.5.

Activity:

Disinfects and cleans in one operation. Bactericidal, fungicidal and virus-inactivating.

Recommended concentration for use in the prevention of hospital cross-infection:

0.75% / onset of action max. 4 hrs.

1.5 % / onset of action max. 1 hr.

Expert Opinion by Prof. Dr. Schubert, Frankfurt 8/87

Prof. Dr. Möse, Graz 9/87

Effective against HBV (hepatitis B) and HIV (AIDS):

1.5% / onset of action max. 2 hrs.

Expert Opinion by Prof. Dr. Frösner, Munich 2/86

Economical combination product. Very good cleaning action. Free of aldehydes, phenols, disinfectant ingredients containing chlorine and phosphates. No unpleasant odor even when using hot water.

Range of application:

Can be used for simultaneous disinfection and cleaning of all washable surfaces and objects throughout the clinical sector and in sanatoria and nursing homes.

Directions for use:

Suitable for any wet-wiping method. Use TANA metering pump for precise dosing. Wipe surface as usual, do not wipe off with water.

Special precautions:

Caustic. Contains quaternary ammonium compounds. Causes acid burns. Must be kept out of the reach of children. Wear protective goggles/face mask. If swallowed, immediately seek medical advice and show the pack or label. In case of contact with the eyes, thoroughly rinse immediately with water and consult a doctor. In case of skin contact, immediately wash off with plenty of water. In accordance with guideline UVV § 7 (Environmental Preventive Regulations), gloves must be worn when handling disinfectants.

This product is intended for industrial use.

Environmental aspects:

This product is subject to the Law on the Environmental Compatibility of Washing and Cleaning Agents and the attached implementing regulations on the biological degradability of raw materials. UBA 0947 0104
